

Supplementary Table 3 Antigen-detecting RDT results with a self-collected NMT swab and a professional-collected NP swab in RT-PCR positive patients, as well as Ct-values, viral load (in descending order), and duration of symptoms per patient. The positive percent agreement between NMT and NP samples on Ag-RDT, and the respective sensitivities compared to RT-PCR are shown.

No.	Ag-RDT		RT-PCR (OP/NP)		Days of Symptoms	Part of self-sampling with relevant deviation
	Self NMT	Prof. NP	Ct value	Viral load ¹		
1	pos. (+++)	pos. (+++)	13.48	9.97 ²	5	rotation
2	pos. (+++)	pos. (+++)	14.74	9.60 ²	1	
3	pos. (+++)	pos. (+++)	14.97	9.53 ²	3	
4	pos. (+++)	pos. (+++)	18.88	9.13 ³	2	
5	pos. (+++)	pos. (+++)	16.41	9.10 ²	2	rubbing
6	pos. (+++)	pos. (+++)	19.67	8.89 ³	2	rotation, rubbing
7	pos. (+++)	pos. (+++)	17.13	8.89 ²	2	rotation, rubbing
8	pos. (+++)	pos. (+++)	17.89	8.66 ²	4	
9	pos. (+++)	pos. (+++)	20.58	8.62 ³	2	rubbing
10	pos. (+++)	pos. (+++)	18.05	8.61 ²	8	
11	pos. (+++)	pos. (+++)	18.08	8.60 ²	2	
12	pos. (+++)	pos. (+++)	18.08	8.60 ²	2	rubbing
13	pos. (+++)	pos. (+++)	20.71	8.59 ³	1	rotation, rubbing
14	pos. (+++)	pos. (+++)	18.80	8.39 ²	4	
15	pos. (+++)	pos. (+++)	21.64	8.31 ³	4	
16	pos. (+++)	pos. (+++)	19.13	8.29 ²	3	rotation, rubbing
17	pos. (+)	pos. (++)	22.74	7.98 ³	4	anterior head tilt
18	pos. (+++)	pos. (+++)	22.85	7.95 ³	1	AN localisation
19	pos. (+++)	pos. (+++)	20.28	7.95 ²	4	rotation, rubbing
20	pos. (+++)	pos. (+++)	22.98	7.91 ³	2	rubbing
21	pos. (++)	pos. (+++)	22.99	7.91 ³	7	
22	pos. (+++)	pos. (+++)	21.83	7.49 ²	2	rotation
23	pos. (++)	pos. (+++)	24.49	7.47 ³	7	unilateral
24	pos. (+)	pos. (++)	21.92	7.46 ²	1	rotation, vertical angle
25	pos. (+++)	pos. (+++)	22.76	7.21 ²	5	rotation
26	pos. (++)	pos. (+++)	26.40	6.90 ³	5	rubbing
27	pos. (+++)	pos. (+++)	26.59	6.85 ³	6	vertical angle
28	pos. (+++)	pos. (+++)	27.13	6.69 ³	-	AN localisation, rubbing
29	neg.	neg.	27.66	6.53 ³	3	rotation
30	pos. (+++)	pos. (+++)	27.82	6.48 ³	1	vertical angle, rotation, rubbing
31	pos. (++)	pos. (+)	25.46	6.41 ²	9	rotation
32	neg.	pos. (++)	25.85	6.29 ²	3	rotation, rubbing
33	pos. (++)	neg.	30.31	5.75 ³	6	
34	neg.	neg.	30.42	4.93 ²	20	rubbing
Sensitivity 31/34 (91.2%)		Sensitivity 31/34 (91.2%)				
Positive percent agreement 96.8%						

¹log₁₀ SARS-CoV2 RNA copies/ml

²TibMolbiol assay, E-gene target

³Roche Cobas SARS-CoV-2 assay (E-gene, T2 target)

Prof., professional sampling; self, self-sampling; AN, anterior nasal; NMT, nasal mid-turbinate; Ag-RDT, Antigen-based rapid diagnostic tool; RT-PCR, Real-time polymerase chain reaction; OP, oropharyngeal; NP, nasopharyngeal; Ct, cycle threshold; No, Number; pos. (+), weak positive; pos. (++) positive; pos. (+++), strong positive; neg., negative