

Supplement 1. Overview of included study patients (n=81) with their study prosthesis (femoral stem, head and cup) grouped by femoral head material. Further columns list additional hip stems (SS=stainless steel), knee and shoulder prostheses implanted either before or within the 10-year study period.

Patient	Stem	Head material	Cup	2nd hip stem (stem material)	1st knee	2nd knee	Shoulder
1	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
2	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	PROFILE (TiAlV)	Genesis I (CoCr/TiAlV)		
3	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE				
4	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE				
5	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	CHARNLEY (SS)	PROFIX (CoCrMo)	PROFIX (CoCrMo)	
6	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
7	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE		PROFIX (CoCrMo)		
8	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
9	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
10	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
11	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE		LCS (CoCrMo)		
12	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE				
13	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE				
14	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	ELITE (CoCrMo)			
15	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
16	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE	LUBINUS SP II (CoCrMo)			
17	Charnley flanged 40	Stainless steel, 316L	Charnley Ogee UHMWPE				
18	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE				
19	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	LUBINUS SP II (CrCo)	LCS (CoCrMo)		
20	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE				
21	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)	PROFIX (CoCrMo)	PROFIX (CoCrMo)	
22	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE				
23	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
24	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
25	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	SPECTRON-EF (CoCrMo)			
26	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	SPECTRON-EF (CoCrMo)			
27	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE				
28	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
29	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE				
30	Spectron EF	CoCrMoMo	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
31	Spectron EF	CoCrMoMo	Reflection XLPE				
32	Spectron EF	CoCrMoMo	Reflection XLPE	LUBINUS SP II (CoCrMo)			
33	Spectron EF	CoCrMoMo	Reflection XLPE				Global Adv. (CoCrMo)
34	Spectron EF	CoCrMoMo	Reflection XLPE	SPECTRON-EF (CoCrMo)			
35	Spectron EF	CoCrMoMo	Reflection XLPE	CHARNLEY (SS)			
36	Spectron EF	CoCrMoMo	Reflection XLPE				
37	Spectron EF	CoCrMoMo	Reflection XLPE	SPECTRON-EF (CoCrMo)			
38	Spectron EF	CoCrMoMo	Reflection XLPE				
39	Spectron EF	CoCrMoMo	Reflection XLPE	LUBINUS SP II (CoCrMo)			
40	Spectron EF	CoCrMoMo	Reflection XLPE	LUBINUS SP II (CrCo)			
41	Spectron EF	CoCrMoMo	Reflection XLPE				
42	Spectron EF	CoCrMoMo	Reflection XLPE	CHARNLEY (SS)			
43	Spectron EF	CoCrMoMo	Reflection XLPE				
44	Spectron EF	CoCrMoMo	Reflection XLPE				
45	Spectron EF	CoCrMoMo	Reflection XLPE	SPECTRON-EF (CoCrMo)			
46	Spectron EF	CoCrMoMo	Reflection XLPE				
47	Spectron EF	CoCrMoMo	Reflection XLPE	EXETER (SS)			
48	Spectron EF	CoCrMoMo	Reflection XLPE	LUBINUS SP II (CoCrMo)	NexGen (TiAlV)	NexGen (TiAlV)	
49	Spectron EF	CoCrMoMo	Reflection XLPE	MS-30 (SS)			
50	Spectron EF	Oxinium	Reflection All-Poly UHMWPE				
51	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	SPECTRON-EF (CoCrMo)	Genesis I (CoCr/TiAlV)		
52	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	PROFILE (TiAlV)			
53	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	Profile/Tri-lock plus (TiAlV)			
54	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	TITAN (TiAlV)			
55	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	SPECTRON-EF (CoCrMo)			
56	Spectron EF	Oxinium	Reflection All-Poly UHMWPE				
57	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	EXETER (SS)			Global Adv. (CoCrMo)
58	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	CHARNLEY (SS)			
59	Spectron EF	Oxinium	Reflection All-Poly UHMWPE				
60	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
61	Spectron EF	Oxinium	Reflection All-Poly UHMWPE				
62	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	SPECTRON-EF (CoCrMo)			
63	Spectron EF	Oxinium	Reflection All-Poly UHMWPE				Delta Xtend (CoCrMo)
64	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	CORAIL (TiAlV)	PROFIX (CoCrMo)	NexGen (TiAlV)	
65	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
66	Spectron EF	Oxinium	Reflection All-Poly UHMWPE	LUBINUS SP II (CoCrMo)			
67	Spectron EF	Oxinium	Reflection XLPE				
68	Spectron EF	Oxinium	Reflection XLPE	EXETER (SS)	NexGen (TiAlV)		
69	Spectron EF	Oxinium	Reflection XLPE	SPECTRON-EF (CoCrMo)	LCS (CoCrMo)		
70	Spectron EF	Oxinium	Reflection XLPE				
71	Spectron EF	Oxinium	Reflection XLPE	SPECTRON-EF (CoCrMo)			
72	Spectron EF	Oxinium	Reflection XLPE	SPECTRON-EF (CoCrMo)			
73	Spectron EF	Oxinium	Reflection XLPE	CHARNLEY (SS)			Global Adv. (CoCrMo)
74	Spectron EF	Oxinium	Reflection XLPE	LUBINUS SP II (CoCrMo)			
75	Spectron EF	Oxinium	Reflection XLPE		PROFIX (CoCrMo)		
76	Spectron EF	Oxinium	Reflection XLPE				
77	Spectron EF	Oxinium	Reflection XLPE	LUBINUS SP II (CoCrMo)			
78	Spectron EF	Oxinium	Reflection XLPE	SPECTRON-EF (CoCrMo)			
79	Spectron EF	Oxinium	Reflection XLPE				
80	Spectron EF	Oxinium	Reflection XLPE	CHARNLEY (SS)			
81	Spectron EF	Oxinium	Reflection XLPE				

Supplement 2. A comparison of the of the blood metal levels ($\mu\text{g/l}$) in patients with only the study prosthesis (pooled) and those with additional prostheses (pooled). An independent samples Mann-Whitney U test was used to compare ranks.

	Groups	N	Medians ($\mu\text{g/l}$)	Mean ranks	p-value
Chromium	Only study prosthesis	26	0.17	40.27	0.85
	Additional prostheses	55	0.20	41.35	
Cobalt	Only study prosthesis	26	0.11	40.13	0.82
	Additional prostheses	55	0.14	41.41	
Zirconium	Only study prosthesis	26	0.07	35.58	0.15
	Additional prostheses	55	0.19	43.56	
Nickel	Only study prosthesis	26	0.17	38.12	0.45
	Additional prostheses	55	0.20	42.36	