

**Supplementary Table 1. List of antibodies for multiplex-immunohistochemistry.** \*CRP stains were done using a Ventana BenchMark Ultra autostainer, all other stains as described in *Methods*

| <i>Antibody</i>              | <i>Clone</i> | <i>Dilution</i> | <i>Product code</i> |
|------------------------------|--------------|-----------------|---------------------|
| <i>Actin (smooth muscle)</i> | 1A4          | 1:500           | DAKO M0851          |
| <i>Caldesmon</i>             | h-CD         | 1:300           | DAKO M3557          |
| <i>CD34</i>                  | QBEnd/10     | 1:50            | DAKO M7165          |
| <i>CD68</i>                  | PG-M1        | 1:100           | DAKO M0876          |
| <i>CK18</i>                  | DC-10        | 1:50            | DAKO M7010          |
| <i>CK7</i>                   | RN7          | 1:200           | NCL-L-CK7-560       |
| <i>CD146</i>                 | UMAB154      | 1:200           | Origene-UM800051    |
| <i>NGFR</i>                  | polyclonal   | 1:500           | Atlas-HPA004765     |
| <i>CRP</i>                   | Y284         | 1:200           | Abcam-ab271830*     |