

Table S1: Detection performance for AIS in a previous validation study.

Statistical Measure	Internal Validation Test	External Validation Test
No. of participants	1050	1072
Sensitivity (%)	92	90
Specificity (%)	88	97
AUROC	0.95	0.97

Table S2: MRI protocols for optional sequences.

Protocol parameters	T1-weighted imaging	T2-weighted imaging
Repetition time (ms)	2090	6310
Echo time (ms)	3	100
No. of excitations	1	2
Matrix size	256 × 240	480 × 336
Flip angle (°)	9	150
Field of view (mm)	210 × 210	210 × 210
Slice thickness (mm)	3.0	4.0
Gap (mm)	4.0	5.0

3T MRI scanner (MAGNETOM Vida, Siemens Healthineers, Erlangen, Germany) with a 64-channel. A head & neck coil was applied for optional sequences.

Table S3: Detection performance for lacunar and non-lacunar AIS.

Statistical Measure	Lacunar AIS	Non-lacunar AIS
No. of participants	63	138
Sensitivity (%)	82.54	92.75
Specificity (%)	89.05	89.05
AUROC (95% CI)	0.91 (0.89-0.93)	0.96 (0.95-0.98)
AUPRC (95% CI)	0.77 (0.65-0.86)	0.91 (0.85-0.95)
Maximum F1 score	0.79	0.87

AIS, acute ischemic stroke; AUROC, area under receiver operating characteristic curve; AUPRC, area under precision-recall curve; CI, confidence interval.