

# New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

## Supplementary Data

**Supplementary Table 1. Digital MIQE2020 guidelines (dMIQE2020) checklist.**

ITEM TO CHECK	PROVIDED	COMMENT
	Y/N	
<b>1. SPECIMEN</b>		
Detailed description of specimen type and numbers	Y	Materials and Methods, references
Sampling procedure (including time to storage)	Y	Materials and Methods, references
Sample aliquotation, storage conditions and duration	Y	Materials and Methods, references
<b>2. NUCLEIC ACID EXTRACTION</b>		
Description of extraction method including amount of sample processed	Y	Materials and Methods, references
Volume of solvent used to elute/resuspend extract	Y	Materials and Methods, references
Number of extraction replicates	Y	Materials and Methods, references
Extraction blanks included?	N	
<b>3. NUCLEIC ACID ASSESSMENT AND STORAGE</b>		
Method to evaluate quality of nucleic acids	Y	Materials and Methods
Method to evaluate quantity of nucleic acid	Y	Materials and Methods
Storage conditions: temperature, concentration, duration, buffer, aliquots	N	
Clear description of dilution steps used to prepare working DNA solution	N	
<b>4. NUCLEIC ACID MODIFICATION</b>		

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

Template modification (digestion, sonication, pre-amplification, bisulfite etc.)	Y	Bisulfite conversion: Materials and Methods, references
Re-extraction performed?	Y	Materials and Methods, references
<b>5. REVERSE TRANSCRIPTION</b>	N	
cDNA priming method + concentration		Not Applicable
One or two step protocol (include reaction details for two step)		Not Applicable
Amount of RNA used per reaction		Not Applicable
Detailed reaction components and conditions		Not Applicable
Estimated copies measured with and without addition of RT		Not Applicable
Manufacturer of reagents used and catalogue number		Not Applicable
Storage of cDNA: temperature, concentration, duration, buffer, aliquots		Not Applicable
<b>6. dPCR OLIGONUCLEOTIDES AND TARGET INFORMATION</b>		
Sequence accession number	N	Commercial kit: Bio-Rad ddPCR MSI RUO Kit
Location of amplicon	N	Commercial kit: Bio-Rad ddPCR MSI RUO Kit
Amplicon length	N	Commercial kit: Bio-Rad ddPCR MSI RUO Kit
Primer and probe sequences (or amplicon context sequence)	N	<p>BAT-25</p> <p>TAAATTAGAACAGGAATTCCAAAGAGACAGCAGTT  GGAACATGAAGAACTAAATTTTCCTCTGCTTTTGGT  TACCACACTTCAAATGACATTCTGCATTTTAACTA  TGGCTCTAAATGCTCTGTTCTCAAAAAAAAAAAAA  AAAAATCAAAAAAAAAACAAAACAAAACTCTTTAGA  GAATCACTCCCACTTACATTCTTGGAGGCGAGGAA  AGCCATGCCCTTTGCCACCTGGTAAGAAAAAGCTCA  GCAAGTCTTCTAAGTCTAGGGCCAACCTCGTCATCCT  CCATGATGGCGGGAGTCACATCTCTTTCTAT</p> <p>BAT-26</p>

# New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

		<p>TTAGAACTCTTATCAGATGATTCCAACCTTTGGACAG  TTTGAAGTACTACTTTTGGACTTCAGCCAGTATATG  AAATTGGATATTGCAGCAGTCAGAGCCCTAACCT  TTTTCAGGTAAAAAAAAAAAAAAAAAGGGTTAA  AAATGTTGAATGGTTAAAAAATGTTTTCATTGACAT  ATACTGAAGAAGCTTATAAAGGAGCTAAAATATTT  TGAAATATTATTATACTTGGATTAGATAACTAGCTT  TAAATGGCTGTATTTTCTCTCCCCTCCTCCACTCCA  CTTTTAACTTTT</p> <p>NR-21</p> <p>TGCTACTCTCTAAAAAAGGCAAGCAGATAAAAGAG  AACACGAAAAATATTCCTACTCCGCATTACACTTT  CTGGTCACTCGCGTTTACAAACAAGAAAAGTGTTG  CTAAAAAAAAAAAAAAAAAGGCCAGGGGAGACATACAT  TTAAATATAAAAATAGAACTGTGCCAGCGACTCCG  GCTGGAATTCTGCTGAAAGGGATGTGTCTTCAGAA  ACCAGGATGTTCTG</p> <p>NR-24</p> <p>GGCTGAGGCAGGAGAATGGCGTGAACCCGGGAGG  CGGACGTTGCAGTGAGCGGAGATTGTGCCATTGCA  TTCCAACCTGGGTGACAGAGTGAGACTCTGTCTCA  CAAAAAAAAAAAAAAAAAATAGGACTGGGCCAGGG  AAGAGAAGAGTTTTTGGAGTCAGGAGGTAAAATTC  AGCAATGGGGATCACGTTAAGGCAGACCTTAAATA  TTCAGATAATTTTCACTAAATAGGGAGCTATTTTA  GGTCTTGAGGAAGGGCAGTGATGGGCAGGAAGGG  GTCATTTG</p> <p>Mono-27:</p> <p>GGCAGGGAAATGGTGGGAACCCAGGGGGTGGAGA  TTGCAGTGAGCTGAGATTGCGCCACTGCACTCCAG  CGTGGGAGACAGAGCAAGACTCTGCCTCAAAAAA  AAAAAAAAAAAAATCCTGGTTTTACTTTTTTTCTTT  TTTAGTTGGCCAAGTGAAATTTGATCCACCCTTAAG  AAAGGAGACAGAACCACATCATGAACCTGTAAGTA  GTATAGCCTTAGAATGTTAGCACTGAAAATAAATG  TTTTAATTTGTTTTTGTGTGATGTATTAGTACTGACC  AAAAAAGCT</p>
Location and identity of any modifications	N	Commercial kit: Bio-Rad ddPCR MSI RUO Kit
Manufacturer of oligonucleotides	N	Commercial kit: Bio-Rad ddPCR MSI RUO Kit
<b>7. dPCR PROTOCOL</b>		
Manufacturer of dPCR instrument and instrument model	Y	Bio-Rad Laboratories QXDx ddPCR System

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

Buffer/kit Catalogue No and manufacturer	Y	Commercial kit: Bio-Rad ddPCR MSI RUO Kit
Primer and probe concentration	N	
Pre-reaction volume and composition (incl. proportion of template added)	Y	6.6ul sample per 22ul reaction
Template treatment (initial heating or chemical denaturation)	N	
Polymerase identity and concentration, Mg <sup>++</sup> and dNTP concentrations	N	
Complete thermocycling parameters	Y	Materials and Methods
<b>8. ASSAY VALIDATION</b>		
Details of optimisation performed		Bio-Rad ddPCR MSI RUO Kit User Guide
Analytical specificity (vs. related sequences)		Bio-Rad ddPCR MSI RUO Kit User Guide
Analytical sensitivity (LoD)	Y	Bio-Rad ddPCR MSI RUO Kit User Guide
Testing for inhibitors (from biological matrix/extraction)	N	Not performed
<b>9. DATA ANALYSIS</b>		
Description of dPCR experimental design	Y	Materials and Methods
Comprehensive details +ve and -ve of controls (whether applied for QC or for estimation of error)	Y	Materials and Methods: positive and negative controls included as part of kit on every plate.
Examples of +ve and -ve experimental results (in supplemental materials)	Y	In results
Description of technical replication	Y	Materials and Methods
Repeatability (intra-experiment variation)	Y	Materials and Methods
Reproducibility (inter-experiment/user/lab etc. variation)	N	

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

Number of partitions measured (average and standard deviation)	Y	Materials and Methods
Partition volume used	N	
Copies per partition ( $\lambda$ or equivalent) (average and standard deviation)	N	
dPCR analysis program (source, version)	Y	QuantaSoft Analysis Pro v1.0 (manual threshold) QX Manager Premium Edition v2.0.0.456 (auto-threshold)
Description of normalisation method	N	Material and Methods: dPCR manual thresholding and auto-thresholding was normalized to the mean cluster densities of the positive controls run on every plate.
Statistical methods used for analysis	N	A Poisson distribution analysis was performed on the dPCR data to calculate copies per ul.
Data submission using RDML	N	Availability of data upon request following all necessary regulatory requirements.

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 2. Validation results obtained for the analysis of FFPE samples using the auto-thresholding algorithm.** Of the over 200 contrived and CRC patient FF, FFPE and cfDNA samples used to develop the auto-thresholding algorithm 60 were FFPE (47 MSI-H and 13 MSS) and 22 were cfDNA (14 MSI-H and 8 MSS), and were used to train the neural network. Another 94 FFPE, 24 MSI-H and 70 MSS, were used to validate the auto-thresholding algorithm. The validation results are shown in the table.

Sample Type	n	Vendor MSI Status*	Concordant Samples by ddPCR MSI	Sensitivity or Specificity	Overall Concordance
FFPE	24	MSI-H	23	95.8%	96.4%
	70**	MSS	67	95.7%	

\*Orthogonal data is immunohistochemistry of mismatch repair proteins

\*\*20 are FFPE of colon tissue from normal patients and presumed MSS

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 3.** Summary of all ddPCR Microsatellite Instability (MSI) RUO kit data that was excluded from the analysis due to a failure or low droplet count of <10,000.

Sample	Sample Type	ddPCR MSI Assay	Reason for Exclusion
Alg-10	FF/FFPE	1	Low droplet count of 4698
Alg-11	FF/FFPE	2	Low droplet count of 8603
Alg-14	FF/FFPE	3	Low droplet count of 2006
Alg-14	FF/FFPE	2	Unsuccessful droplet generation, no droplets detected
Alg-16	FF/FFPE	1	Low droplet count of 1396
Alg-16	FF/FFPE	2	Unsuccessful droplet generation
Alg-23	FF/FFPE	1	Low droplet count of 3699
Alg-28	FF/FFPE	2	Unsuccessful droplet generation, no droplets detected
Alg-34	FF/FFPE	3	Low droplet count of 3342
Alg-51	FF/FFPE	2	Unsuccessful droplet generation
Alg-58	FF/FFPE	2	Low droplet count of 3047
Alg-58	FF/FFPE	3	Low droplet count of 5056
Alg-61	FF/FFPE	3	Low droplet count of 2816
Alg-67	FF/FFPE	2	Low droplet count of 5493
Alg-70	FF/FFPE	3	Low droplet count of 5938
Alg-73	FF/FFPE	1	Insufficient plate seal before thermal cycling
Alg-91	FF/FFPE	3	Low droplet count of 7958
Alg-93	FF/FFPE	1	Low droplet count of 675
Ras-95	cfDNA	3	Cluster mirroring

# New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 4.** Concordance of FF and FFPE tissues from colorectal cancer patients between auto-thresholded ddPCR MSI with the Promega MSI Analysis System v1.2 (A). There were twelve samples that were not concordant for all markers: Samples <sup>†</sup>Alg-96; <sup>\*</sup>Alg-25, Alg-67; <sup>‡</sup>Alg-99, <sup>†</sup>Alg-45, <sup>¶</sup>Alg-30, Alg-43, Alg-53, Alg-63; <sup>\*\*</sup>Alg-17; <sup>††</sup>Alg-88; <sup>§</sup>Alg-28. For Alg-28, the same 4 markers were not unstable for both ddPCR MSI and the MSI Analysis System v1.2. Eleven samples were re-tested with the MSI Analysis System v1.2 (Alg-67 not re-tested) (B). Seven of these were concordant compared to ddPCR MSI, while four remained discordant by marker: Sample <sup>\*</sup>Alg-25 & Alg-67, <sup>‡</sup>Alg-99, <sup>¶</sup>Alg-30, <sup>§</sup>Alg-28.

**A**

		Promega MSI Analysis (tissue DNA), Historical					
		MSS		MSI-H			
		0	1	2	3	4	5
ddPCR MSI Auto-threshold (tissue DNA)	MSS	0	83				
		1	1 <sup>†</sup>				
		2	2 <sup>*</sup>				
	MSI-H	3	1 <sup>‡</sup>				
		4				1 <sup>§</sup>	
		5	1 <sup>†</sup>	4 <sup>¶</sup>	1 <sup>*</sup>	1 <sup>†</sup>	7

**B**

		Promega MSI Analysis + Discordant Repeat Result (tissue DNA), Historical					
		MSS		MSI-H			
		0	1	2	3	4	5
ddPCR MSI Auto-threshold (tissue DNA)	MSS	0	83				
		1		1			
		2	2 <sup>*</sup>	1 <sup>‡</sup>			
	MSI-H	3					
		4					1 <sup>§</sup>
		5			1 <sup>¶</sup>		13



# New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 5.** Concordance of FF and FFPE tissues from colorectal cancer patients between manual thresholded ddPCR MSI with the Promega MSI Analysis System v1.2 (A). There were fourteen samples that were not concordant for all markers: Sample <sup>†</sup>Alg-1, Alg-67, Alg-96, Alg-99; \*Alg-25; <sup>‡</sup>Alg-45; <sup>¶</sup>Alg-30, Alg-43; \*\*Alg-53, Alg-63; <sup>‡</sup>Alg-17; <sup>§</sup>Alg-28; <sup>††</sup>Alg-88; <sup>#</sup>Alg-75. For Sample Alg-28, the same 4 markers were not unstable for both ddPCR MSI and the MSI Analysis System v1.2. Twelve of the fourteen samples were re-tested with the MSI Analysis System v1.2 (B). Sample Alg-1 and Alg-67 not retested. Seven of these were concordant compared to ddPCR MSI, while five remained discordant by marker: Sample <sup>†</sup>Alg-1, Alg-67; \*Alg-25; <sup>¶</sup>Alg-30, <sup>‡</sup>Alg-75; <sup>§</sup>Alg-28, Alg-43. The concordance by MSI status was high at 101 of 102 (C).

**A**

		Promega MSI Analysis (tissue DNA), Historical					
		MSS		MSI-H			
		0	1	2	3	4	5
ddPCR MSI Manual Threshold (tissue DNA)	MSS	0	82				
		1	4 <sup>†</sup>				
		2	1 <sup>*</sup>				
	MSI-H	3					1 <sup>#</sup>
		4		2 <sup>¶</sup>		1 <sup>§</sup>	
		5	1 <sup>‡</sup>	2 <sup>*</sup>	1 <sup>‡</sup>	1 <sup>†</sup>	6

**B**

		Promega MSI Analysis + Discordant Repeat Result (tissue DNA), Historical					
		MSS		MSI-H			
		0	1	2	3	4	5
ddPCR MSI Manual Threshold (tissue DNA)	MSS	0	82				
		1	2 <sup>†</sup>	2			
		2	1 <sup>*</sup>				
	MSI-H	3					1 <sup>‡</sup>
		4			1 <sup>¶</sup>		2 <sup>§</sup>
		5					11

**New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability**

**C**

		<b>Promega MSI Analysis system + Discordant Repeat Test</b>	
		<b>MSS</b>	<b>MSI-H</b>
<b>ddPCR MSI Manual Threshold</b>	<b>MSS</b>	86	0
	<b>MSI-H</b>	1	15

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 6.** Comparison of ddPCR MSI with auto-thresholding for 14 FF and FFPE samples when tested at an input per well of 18, 5 and 1 ng/well, measured by Qubit dsDNA BR Assay. Seven (7) of the samples were determined as MSS and 7 as MSI-H by the Promega MSI Analysis System. The ddPCR MSI results reflect data for replicate 2 of a total of 2 replicates.

Sample	Promega MSI Analysis					ddPCR MSI - 18 ng/rxn					ddPCR MSI - 5 ng/rxn					ddPCR MSI - 1 ng/rxn					Status			
	RAT-25	RAT-26	NR-21	NR-24	Mono-27	Status	BAT-25	BAT-26	NR-21	NR-24	Mono-27	Status	RAT-25	RAT-26	NR-21	NR-24	Mono-27	Status	BAT-25	BAT-26		NR-21	NR-24	Mono-27
Alg-1	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	+	-	MSS	-	-	-	-	-	MSS
Alg-2	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	+	MSS
Alg-14	-	-	-	-	-	MSS	-	-	na	na	na	na	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-16	-	-	-	-	-	MSS	na	na	na	na	-	na	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-67	-	-	-	-	-	MSS	-	+	na	na	-	na	-	+	-	-	-	MSS	-	-	-	-	+	MSS
Alg-87	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-100	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-28	+	+	+	+	+	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	+	MSI-H
Alg-30	+	+	+	-	-	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-33	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-43	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-58	+	+	+	+	+	MSI-H	+	+	na	na	na	na	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-83	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-89	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H

-	stable
+	unstable
na	failed

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 7.** Comparison of ddPCR MSI with manual thresholding for 14 FF and FFPE samples when tested at an input per well of 18, 5 and 1 ng/well, measured by Qubit dsDNA BR Assay. Seven (7) of the samples were determined as MSS and 7 as MSI-H by the Promega MSI Analysis System. The ddPCR MSI results reflect data for replicate 1 of a total of 2 replicates.

		Promega MSI Analysis					ddPCR MSI - 18 ng/rxn					ddPCR MSI - 5 ng/rxn					ddPCR MSI - 1 ng/rxn								
Sample		BAT-25	BAT-26	NR-21	NR-24	Mono-27	Status	BAT-25	BAT-26	NR-21	NR-24	Mono-27	Status	BAT-25	BAT-26	NR-21	NR-24	Mono-27	Status	BAT-25	BAT-26	NR-21	NR-24	Mono-27	Status
Alg-1		-	-	-	-	-	MSS	-	+	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-2		-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-14		-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-16		-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-67		-	-	-	-	-	MSS	-	+	-	-	-	MSI-H	-	+	-	-	-	MSS	-	-	-	-	-	MSS
Alg-87		-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-100		-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-28		+	+	+	+	+	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H	-	+	+	+	-	MSI-H
Alg-30		+	+	+	-	-	MSI-H	+	+	+	-	+	MSI-H	+	+	-	-	+	MSI-H	+	+	+	+	-	MSI-H
Alg-33		+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-43		+	+	+	+	+	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H
Alg-58		+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-83		+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-89		+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H

-	stable
+	unstable
na	failed

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 8.** Comparison of ddPCR MSI with manual thresholding for 14 FF and FFPE samples when tested at an input per well of 18, 5 and 1 ng/well, measured by Qubit dsDNA BR Assay. Seven (7) of the samples were determined as MSS and 7 as MSI-H by the Promega MSI Analysis System. The ddPCR MSI results reflect data for replicate 2 of a total of 2 replicates.

Sample	Promega MSI Analysis					ddPCR MSI - 18 ng/rxn					ddPCR MSI - 5 ng/rxn					ddPCR MSI - 1 ng/rxn								
	RAT-25	RAT-26	NR-21	NR-24	Mono-27	Status	BAT-25	BAT-26	NR-21	NR-24	Mono-27	Status	RAT-25	RAT-26	NR-21	NR-24	Mono-27	Status	RAT-25	RAT-26	NR-21	NR-24	Mono-27	Status
Alg-1	-	-	-	-	-	MSS	-	+	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-2	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-14	-	-	-	-	-	MSS	-	-	na	na	na	na	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-16	-	-	-	-	-	MSS	na	na	na	na	-	na	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-67	-	-	-	-	-	MSS	-	+	na	na	-	na	-	+	-	-	-	MSS	-	-	-	-	-	MSS
Alg-87	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-100	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS	-	-	-	-	-	MSS
Alg-28	+	+	+	+	+	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H
Alg-30	+	+	+	-	-	MSI-H	+	+	+	-	+	MSI-H	+	+	-	+	+	MSI-H	+	+	-	-	+	MSI-H
Alg-33	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-43	+	+	+	+	+	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H	+	+	+	+	-	MSI-H
Alg-58	+	+	+	+	+	MSI-H	+	+	na	na	na	na	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-83	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H
Alg-89	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H	+	+	+	+	+	MSI-H

-	stable
+	unstable
na	failed

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 9. Intra- and inter-run reproducibility of sample testing.** Nine runs for cohort 1 (ALGECOLS) were analyzed for intra-and inter-run reproducibility by analyzing the positive controls run in duplicate for each run. Intra-run analysis of the positive control resulted in a mean difference in mutant fractional abundance ranging from 5.8% to 13.2% across the 5 target loci. The mean difference in wild type concentration for the 3 assays ranged from 2.8% to 3.5%. Inter-run analysis of the positive control resulted in a mutant fractional abundance CV ranging from 4.0% to 9.6% (n = 18) across the 5 target loci. The CV of the wild type concentration ranged from 2.0% to 2.5% (n = 18).

<b>Marker</b>	<b>Mean Difference in Mutant Fractional Abundance between Duplicates</b>	<b>Mean Mutant Fractional Abundance (n = 18)</b>	<b>CV of Mutant Fractional Abundance (n = 18)</b>
BAT-25	8.1%	7.1%	5.8%
BAT-26	11.4%	6.9%	7.8%
NR-21	13.2%	6.9%	9.6%
NR-24	5.8%	7.5%	4.0%
Mono-27	10.6%	6.0%	7.9%
<b>Wild Type</b>	<b>Mean Difference in Wild Type Concentration between Duplicates</b>	<b>Mean Wild Type Concentration (copies/<math>\mu</math>L) (n = 18)</b>	<b>CV of Wild Type Concentration (n = 18)</b>
Assay 1	3.4%	170.0	2.4%
Assay 2	3.5%	169.5	2.5%
Assay 3	2.8%	172.0	2.0%

## New Multiplex Droplet Digital PCR for highly sensitive detection of microsatellite instability

**Supplementary Table 10. Comparison of duplicate analysis.** For cohort 1, 102 samples were run in duplicate within the same run. Eighty-eight (88) of the 102 had valid results for both replicates. Among these 88 samples, the mean difference in wild type concentration for the 3 assays ranged from 4.3% to 4.5%. The marker results for these 88 samples had a 100% match between replicates for 86 of the samples. For the remaining 2 samples, detection of one marker changed from positive to negative because the mutant fractional abundance of these markers was near the limit of detection. These discordant marker results did not change the MSI status determination for these samples. Of the 88 samples, the ddPCR MSI assay determined 16 were MSI-H. The mean difference in mutant fractional abundance for these 16 samples ranged between 2.1% and 6.5% for the 5 target loci.

	Mean Difference between Wild Type Concentration of Duplicates			Standard Deviation of Difference between Wild Type Concentration of Duplicates		
Wild Type	MSS (n = 72)	MSI-H (n = 16)	All Samples (n = 88)	MSS (n = 72)	MSI-H (n = 16)	All Samples (n = 88)
Assay 1	4.2%	4.7%	4.3%	5.3%	6.2%	5.4%
Assay 2	4.5%	4.4%	4.5%	6.0%	5.6%	5.9%
Assay 3	4.4%	3.8%	4.3%	5.7%	4.9%	5.5%
Marker	Mean Difference in Mutant Fractional Abundance of MSI-H  n = 16			Standard Deviation of Difference in Mutant Fractional Abundance of MSI-H  n = 16		
BAT-25	5.5%			10.0%		
BAT-26	4.1%			7.3%		
NR-21	2.1%			3.1%		
NR-24	6.5%			12.8%		
Mono-27	5.8%			11.5%		