Supplementary Note: Comparison of standard and ultra-rapid ddPCR thermal conductance and surface area to volume ratios

Here, we estimate the thermal conductances of the ddPCR 96-well plate well used in standard ddPCR and the stainless-steel capillary used in ultra-rapid ddPCR. Thermal conductance, C, is the rate at which a material of a given thickness conducts heat. It is calculated as C = K/T and measured in Watts/(meter^{2*}Kelvin), i.e., W/(m^{2*}K), where K is the thermal conductivity (a heat transfer constant specific to a given material) measured in Watts/(meter*Kelvin) and T is the thickness of the material, measured in meters.

The thermal conductivity of polypropylene, the plastic used in the ddPCR 96well plates, has been measured to be 0.22 W/(m*K) ¹, and the thermal conductivity of stainless-steel 304, the material of our capillaries, is approximately 15 W/(m*K) ². Since the thickness of the ddPCR plate wells was not available to us, we measured the thickness of 33 wells of a ddPCR 96-well plate at their thinnest point using a micrometer and obtained an estimate of 0.30 \pm 0.07 mm. The thickness of the stainless-steel capillary is 0.10 \pm 0.04 mm per the manufacturer. Thus, the thermal conductance of the ddPCR plate well is approximately 7.3 ×10² W/(m²*K) and the thermal conductance for the stainless-steel capillary is approximately 1.5 ×10⁵ W/(m²*K), which is an ~ 204-fold difference.

We also estimated the surface area to volume ratio of a standard ddPCR reaction occurring in a ddPCR 96-well plate and an ultra-rapid ddPCR reaction occurring in a stainless steel capillary. We estimated the volume of the slanted portion of an individual well of a ddPCR 96-well plate as a truncated cone using bottom and top diameters of 2.16 mm and 5.60 mm, respectively, and a height of 11.10 mm as estimated from the manufacturer's blueprint, which yielded a volume of 139.84 μ L. Using these values, we also calculated the surface area of the slanted portion of the ddPCR plate well to be 136.92 mm². Thus, the surface area to volume ratio of a standard ddPCR well is 0.98 mm²/ μ L. For the stainless-steel capillary, we calculated its volume as a cylinder with length 57.15 mm and an inner diameter of 1.067 mm (per the manufacturer), which yielded 51.10 μ L. Using these values, we also calculated the surface area contacting the heated water during ultra-rapid ddPCR, to be 191.57 mm². Thus, the surface area to volume ratio of the capillary is 3.75 mm²/ μ L, which is 3.8-fold higher than that of the ddPCR plate well.

References

- 1 Antonella, P. & Domenico, A. in *Polypropylene* (eds Wang Weiyu & Zeng Yiming) Ch. 3 (IntechOpen, 2019).
- 2 Graves, R. S., Kollie, T. G., McElroy, D. L. & Gilchrist, K. E. The thermal conductivity of AISI 304L stainless steel. *International Journal of Thermophysics* **12**, 409-415 (1991).