

Supplementary Table 1. Reference laboratory blinded SARS-CoV-2 quality assessment panel.

Gamma irradiated culture supernatant of SARS-CoV-2: Dilution	Reference laboratory RdRP gene RT-PCR: Cycle threshold	EZ1 Advanced extraction/Coronavirus Typing assay SARS-CoV-2: Semi-quantitative take-off value	MT-PREP extraction/Coronavirus Typing assay SARS-CoV-2: Semi-quantitative take-off value
1:1000	29	9.48	9.41
1:1000	29	10	9.57
1:4000	31	12.39	12
1:4000	31	12.86	11.98
1:16000	33	14.76	14.33
1:16000	33	14.43	14.65
1:64000	35	17.5	16.2
1:64000	35	16.9	16.39
1:128000	36	20.11	17.95
1:128000	36	20.19	17.98
1:256000	37	19.25	19.47
1:256000	37	18.6	18.56

Supplementary Table 2. Other respiratory pathogens identified as part of routine clinical care from samples collected for SARS-CoV-2 testing during study period.

Respiratory pathogen	Detected	Total tested
SARS-CoV-2	54 (2.0%)	2673
HCoV-229E	13 (0.4%)	2673
HCoV-HKU1	26 (0.9%)	2673
HCoV-NL63	31 (1.2%)	2673
HCoV-OC43	49 (1.8%)	2673
Influenza A	25 (2.7%)	914
Influenza B	2 (0.2%)	914
RSV	16 (1.8%)	914
Parainfluenza 1	1 (1.6%)	61
Parainfluenza 2	0 (0%)	61
Parainfluenza 3	1 (1.6%)	61
Parainfluenza 4	0 (0%)	61
Rhinovirus/enterovirus	13 (13.5%)	96
Adenovirus	4 (4.2%)	96
HMPV	8 (22.9%)	35
<i>Mycoplasma pneumoniae</i>	2 (2.1%)	96

Abbreviations: CoV, coronavirus; HCoV, human coronavirus; hMPV, Human Metapneumovirus; MERS, Middle Eastern respiratory syndrome; Orf, open reading frame; RSV, Respiratory syncytial virus; SARS, Severe acute respiratory syndrome.