

Supplementary data 1: Multi-echo 3D FLASH acquisition parameters for Siemens Magnetom PRISMA MRI

	PDw	T1w	MTw
TR [ms]	24.5	24.5	24.5
Flip angle [°]	6	21	6
Bipolar gradient echoes/TE [ms]	8/TE 2.34 – 18.72	8/TE 2.34 – 18.72	6/TE 2.34 – 14.04
Off-resonance Gaussian MT pulse	N/A	N/A	FA: 220° Frequency offset: 2 [kHz]
Bandwidth [Hz/Px]	465	465	465