

Supplemental Methods

Hydrophobicity index

To calculate the total plasma bile acid hydrophobicity index (HI) we used the previously published HI for each individual bile acid¹ and multiplied them by the relative concentration of each bile acid.

Formula:

Total plasma bile acid HI=

((bile acid 1 concentration * bile acid 1 HI)/ total bile acid amount)

+ (bile acid 2 concentration * bile acid 2 HI)/ total bile acid amount)

+ ...

+ (bile acid 14 concentration * bile acid 14 HI)/ total bile acid amount))

Reference

1. Heuman DM. Quantitative estimation of the hydrophilic-hydrophobic balance of mixed bile salt solutions. *J Lipid Res.* 1989;30(5):719-730. doi:[https://doi.org/10.1016/S0022-2275\(20\)38331-0](https://doi.org/10.1016/S0022-2275(20)38331-0)